


USING ADDITIONAL FLEXTEND® GLOVE RINGS

GLOVE PROFILE: FRONT

The additional ring located on the front of the FLEXTEND® glove, circled in **BLUE**, is for use with **both** the FLEXTEND® and FLEXTEND®-AC Kit.

- **FIG. 1** - FLEXTEND® USERS: Placing the blue or black elastic resistance band through this ring will provide increased thumb resistance.
- **FIG. 2** - FLEXTEND®-AC Kit: Used to attach AC-Kit resistance bands. (See FLEXTEND® AC-Kit Manual)



GLOVE PROFILE: BACK & SIDES

The additional ring located on the back and sides of the glove marked in **BLUE** and **RED** are for use with the FLEXTEND® AC-Kit.

- **FIG. 3** - FLEXTEND®-AC Kit: Used to attach AC-Kit resistance bands. (See FLEXTEND® AC-Kit Manual)



NOTE: If you purchased the AC-KIT Upper Extremity Exercise Attachment System, you should have received an Instruction Manual in the Kit itself. If you have questions or would like to upgrade your FLEXTEND® purchase by adding the AC-Kit, please contact original place of purchase.